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Your opinion matters! Visit LearningResources.com to write a product review or to find a store near you.

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ACTIVITIES

Order Up! (K–2)

Fill equal-size containers with different amounts of colored water. Then, line up the containers based on the capacity of water the containers hold.

Estimate Equivalents (2+)

Explore a variety of the included measurement tools. Then, estimate how many of the smaller containers will fit into the larger containers. Measure the actual amounts by filling the containers with water. Record how many units were needed to fill each container (e.g., How many quarts are in a gallon?). Extension activity: Use a variety of everyday containers (e.g., soda bottle, soup can, milk jug) and find their equivalents.

Cooking with Equivalents (2+)

Explore a variety of recipes. Make one of the recipes in class. Try doubling the recipe or making half the servings. For more advanced students, have them convert the recipe from standard to metric and vice versa.

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Your opinion matters! Visit
LearningResources.com
to write a product review or to
find a store near you.

Like
us on

© Learning Resources, Inc., Vernon Hills, IL, US
Learning Resources Ltd., Bergen Way,
King's Lynn, Norfolk, PE30 2JG, UK
Please retain the package for future reference.
Made in China. LRM2420-GUD

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
3+ | 8+

ages
años • ans • jahre

Precision School Balance with Weights

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Guía de actividades • Guide d'activités • Spielvorschläge

CAUTION:

Federal law requires us to advise that THE BRASS WEIGHTS in this educational product may contain lead which may be harmful if swallowed. We stand behind the safety of all of our products. For further information please contact us at QA@learningresources.com or 800-222-3909.

WARNING:
CHOKING HAZARD - Small parts.
Not for children under 3 years.

Help students explore, develop, and extend their understanding of mathematical and scientific concepts including inequalities, mass, comparisons, the four operations, and order. The Precision School Balance comes with a 10-piece gram set including two 1-gram masses, two 2-gram masses, two 5-gram masses, two 10-gram masses, one 20-gram mass, one 50-gram mass.

Balance Components

- The damping button helps make quick and accurate measurements. The button steadies the swing arm and prevents the pans from rocking during a measurement.
- The calibration slider allows you to zero the balance.
- Marked calibrations demonstrate when the balance is level and zeroed.
- The balance contains a case with a lid to hold the masses. Storing your masses in this compartment will reduce the likelihood of misplacing them.

Use the balance in directed activities such as those listed below or simply place it in a learning center for students to explore. Either way, this balance will expand and enrich the learning opportunities in your classroom.

Exploring the Precision School Balance

Collect a variety of classroom objects such as pencils, erasers, paper clips, scissors, counters, or blocks. Encourage students to place objects on either side of the balance and observe the results. Ask questions that will increase your students' explorations.

You might ask:

- Are there any objects that balance each other?
- Does size always tell which object has the greatest mass?
- Can a small object have greater mass than a large object?

After students have explored with objects on either side of the balance, encourage them to explore measurement with the unit masses provided with the balance. Students can place one object on one side of the balance and use masses to find out if an object has greater, lesser, or close to equal the mass in a unit weight.

Equal Masses

Once students have explored measurements that are greater, lesser, or close to the actual mass, encourage them to find the actual mass. Introduce them to the calibration lines that show when the balance is level. Ask, "Which line should the pointer be on for objects to be balanced?" (Answer: The center line.)

Discuss how to determine an object's mass. Provide examples and allow students to discover the most efficient method. You might prompt their explorations with questions such as:

- Is it best to start with the small unit masses to measure an object?
- Is it best to start with the large unit masses to measure an object?

Ordering Masses

Encourage students to arrange, by mass, the items they have collected. Students should choose one item, determine its mass, and record its mass in a data table. After doing this for all the objects, students should arrange them in order from the least mass to the greatest mass.

Mass Hunt

Play a game of Mass Hunt! Hold up an object and ask students to identify it, and then determine its mass. Challenge them to find another object in the classroom that has a similar mass. Give each student the opportunity to choose and reject items until they agree on one item that they think has the mass that most closely equals that of the object you selected. Discuss students' choices and then determine the masses of the most popular choices to see which one is the best answer.

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